

Spine Lumbar Torso with Pelvis, Muscles and Nerves in Soft Tissue and Skin

The Spine Lumbar Torsos with pelvis, muscles and nerves in soft tissue with skin were specially developed for surgeons to practice various MIS-procedures on the dorsal access to the spine. The Torsos can be used for screwing and plating, placing rods and external fixations as well as for imaging. Available in various sizes, materials and pathologies.

Features and Benefits

- Soft tissue firmly surrounds spine model
- Red soft tissue for improved contrast
- Near to human-like feeling
- Easy to cut the skin
- Cortical shell and cancellous bone structure
- Radio-opaque coating for CT-imaging
- Fulfils all anatomical characteristics
- Rigid hold of spine, allowing corrective procedures



Spine Lumbar Torso with Pelvis, Muscles and Nerves in Soft Tissue and Skin

Available sizes from L4 to L5 intact
from L1 to L5 intact
from T10 to L5 with or without scoliosis

Materials Spine: LD (Solid Foam Density) or LSS (Cancellous Bone Density)
Disks: Flexible Foam
Muscles: Flexible Foam
Nerves: Flexible Foam
Skin: PU Material

Art Numbers **Models in solid foam density**
LD9311.10.FI L4-L5 intact
LD9312.FI L1-L5 intact
LD9313.FI T10-L5 intact
LD9361.FI T10-L5 with scoliosis

Models in cancellous bone density
LSS9311.10.FI L4-L5 intact
LSS9312.FI L1-L5 intact
LSS9313.FI T10-L5 intact
LSS9361.FI T10-L5 with scoliosis

contact:
sales@synbone.com

To obtain best training results we recommend using SYNBONE training models.
Contact us for our complete SPINE product portfolio.

SYNBONE AG • Tardisstrasse 199 • 7205 Zizers • Switzerland
Phone +41 81 300 02 80 • sales@synbone.com • www.synbone.com

SYNBONE SDN BHD • No. 40-41, Jalan i-Park 1/3 • 81000 Bandar • Indahpura • Kulai, Johor • Malaysia
Phone +60 7660 8220 • sales@synbone.com • www.synbone.com